

Handheld Dry Chemistry Analyzer

DCA-25

HEALTH MANAGER AT YOUR FINGERTIPS



RAPID

Test result in **60 s**



COMPACT AND PORTABLE

151.5*81*28 mm 230 g



INTELLIGENT

Available for cloud health app

Functions



Key Parameters

Methodological	Dry Chemistry
Detection Time	60s
Sample Type	Whole blood, serum, plasma, peripheral blood, etc.
Weight	230 g
Size	151.5mm (length) × 81mm (width) × 28 mm (height)
Data Transmission	Bluetooth, USB

Test Items

Creatinine

Matrix Labs Diagnostics Private Ltd.,

S.No-114/13A, 13B, 12B, 12A2, Ground Floor,
Ponnamman Koil Street, Noombal,
Chennai-600077, Thiruvallur District, Tamilnadu, India.
+91-044-26793885, Toll free : 18001201713,
info@matrixlabs.in, www.matrixlabs.in

First Floor, New Anna, Chennai
600035 P.H. O.

